

Safety Injection Series Solutions

S ETE

Better care, Better life.

LIIIIIIIIIIIIIII



Sterile Syringe For Single Use

Features:

- Central nozzle type or eccentric nozzle type
- Luer-slip type or luer-lock type
- Two-piece type or three-piece type
- Soft medium container or hard medium container
- With needle or without needle.

Specifications: 1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50ml, 60ml

Hypodermic Needle

Features:

- The needle made of high quality austenitic stainless steel. adopts the international popular thin-walled tube design, with large inner diameter and high flow rate.
- Fast needle entry, less pain and less tissue damage;
- Needle holder with different color to identify specifications.
- Specifications: 0.45 mm to 1.8mm





Disposable Insulin Syringe

Features:

- Ultra-micro needle tube, minimizing the injection pain.
- Barrel and plunger: Polypropylene
- Piston: Natural Latex/Latex Free

Specifications: 1ml

Sterile Syringes for Fix-dose Immunization for Single Use

Features:

- The biggest feature is The Fixed Dose
- Medical grade high transparent PP
- Clear barrel easy to observe liquid.the ink of scale has strong adhesion
- High quality rubber with excellent sealing and sliding properties.

Specifications: 0.1ml, 0.5ml, 1ml, 2ml, 2.5ml, 3ml



FFEE

Auto-Disable Syringe

Features:

- The self-destruction function will be automatically started after injection, effectively preventing secondary use
- ◆ Special design to ensure the product micro-residue amount, and the drug solution dosage is effectively used. The national standard is ≤0.07ml

Specifications: 0.1ml, 0.5ml, 1ml, 2ml, 2.5ml, 3ml, 5ml

Auto-Retractable Syringe

Features:

- The injection needle is completely pulled back into the sheath to prevent the risk of needle sticks
- The special structure design enables the conical connector to drive the injection needle assembly to completely retract into the sheath, effectively preventing the risk of needle sticks for medical staff.
- Dynamic interference fit between components to ensure that the product does not leak.



Specifications: 1ml, 3ml, 5ml, 10ml

Needle-Free Injection System

Features:

- Painless injection to relieve the psychological pressure of patients
- Subcutaneous diffusion technology to improve drug absorption rate
- Needle-free injection to avoid needle stick injuries of medical staff
 - Protect the environment and solve the medical waste recycling problem of traditional injection devices





Needle injection

Needle free injection

inning and announce





江西三鑫医疗科技股份有限公司

JIANGXI SANXIN MEDTEC CO.,LTD. 地址: 江西南昌小蓝经济技术开发区富山大道999号 Add: No.999 Fushan Road, Xiaolan Economic Development Zone 330200 Nanchang,Jiangxi P.R. China 电话 (Tel) : 86-791-85950275 Fax: 86-791-85988819 邮箱 (E-mail) : sales@sanxin-med.com 网址 (Web) : www.sanxin-med.com





Sterile Safety Needle for Single Use

江西三鑫医疗科技股份有限公司

JIANGXI SANXIN MEDTEC CO.,LTD. 地址: 江西南昌小蓝经济技术开发区富山大道999号 Add: No.999 Fushan Road, Xiaolan Economic Development Zone 330200 Nanchang,Jiangxi P.R. China 电话 (Tel) : 86-791-85950275 Fax: 86-791-85988819 邮箱 (E-mail) : sales@sanxin-med.com 网址 (Web) : www.sanxin-med.com



FEATURES:

- Safety Mechanism activated by one hand.
- Audible click confirms completion of safety device.
- No change in puncture technique.
- Ultra sharp three-beveled needle.
- Siliconization of needle minimizes pain from patient.
- Colored cap conform to ISO6009.
- Medical Grade Material with latex free, DEHP free.
- Compatible for both Lure slip and lure lock syringe.

PRODUCT RANGE OVERVIEW:



| Product | Size | Lure Lock (LL) /Lure Slip (LS) | Needle Gauge | Needle Length |
|----------------------------|---------|-----------------------------------|--------------|---------------|
| Syringe with Safety Needle | 1ML | LL/LS | 18G-27G | 16mm-38mm |
| Syringe with Safety Needle | 2/2.5ML | LL/LS | 18G-27G | 16mm-38mm |
| Syringe with Safety Needle | 3ML | LL/LS | 18G-27G | 16mm-38mm |
| Syringe with Safety Needle | 5ML | LL/LS | 18G-27G | 16mm-38mm |
| Syringe with Safety Needle | 10ML | LL/LS | 18G-27G | 16mm-38mm |
| Syringe with Safety Needle | 20ML | LL/LS | 18G-27G | 16mm-38mm |
| Syringe with Safety Needle | 30ML | LL/LS | 18G-27G | 16mm-38mm |
| Syringe with Safety Needle | 50ML | LL/LS | 18G-27G | 16mm-38mm |
| Syringe with Safety Needle | 60ML | LL/LS | 18G-27G | 16mm-38mm |

SAFETY NEEDLE INSTRUCTION:



