



-B4-



-B1-

**江西三鑫医疗科技股份有限公司**

JIANGXI SANXIN MEDTEC CO.,LTD.

地址: 江西南昌小蓝经济技术开发区富山大道999号

Add: No.999 Fushan Road, Xiaolan Economic Development

Zone 330200 Nanchang, Jiangxi P.R. China

电话 (Tel) : 86-791-85950275 Fax: 86-791-85988819

邮箱 (E-mail) : sales@sanxin-med.com

网址 (Web) : www.sanxin-med.com




**Sansin**



**Dialysis Catheter Needle  
for Single Use**

- BPA Free
- DEHP Free
- Material FEP




Y  
Type



Straight  
Type

## Features



## Advantages

- Effectively decrease the vascular access damage compare to AVF.
- Greatly reduce the rate of hematoma and errhysis.
- Safety Needle tip device help to prevent the injury to nurse and doctors.
- Flexible FEP catheter improve the compatibility.
- Easy to dismantle.

## Better choice for patients

- Poor fistula condition patients.
- Allergic to metal needle patients.
- Night-time dialysis patients.
- Patients during Fistula maturation period.
- Patients with anxiety and irritability (dementia patient)
- Patients with infectious diseases (HBV/HCV/HIV)



## Specifications and models

Model	Structure		Safety Device	Sleeve	Wing	Clamp	Cap
	Straight Type	Y Type					
B1	✓						
B2	✓				✓	✓	
B3	✓				✓	✓	
B4	✓				✓		
B5	✓		✓				
B6	✓		✓		✓	✓	
B7	✓		✓		✓	✓	
B8	✓		✓		✓		
Y1		✓					✓
Y2		✓				✓	✓
Y3		✓				✓	✓
Y4		✓					
Y5		✓				✓	
Y6		✓				✓	
Y7		✓	✓	✓			✓
Y8		✓	✓	✓		✓	✓
Y9		✓	✓	✓		✓	✓
Y10		✓	✓	✓			
Y11		✓	✓	✓		✓	
Y12		✓	✓	✓		✓	