Sansin



江西三鑫医疗科技股份有限公司

JIANGXI SANXIN MEDTEC CO.,LTD.

地址: 江西南昌小蓝经济技术开发区富山大道999号

Add: No.999 Fushan Road, Xiaolan Economic Development

Zone 330200 Nanchang, Jiangxi P.R. China

电话 (Tel): 86-791-85950275 Fax: 86-791-85988819

邮箱 (E-mail) : sales@sanxin-med.com 网址 (Web) : www.sanxin-med.com





1. Equipments

- 1-1. RO Water Room
- 1-2. CDS/Concentration Room
- 1-3. Hemodialysis Machine
- 1-4. Dialysis Powder Mixer

2. Hemodialysis Consumable

- 2-1. Wetted Consumable Room
 - 2-1-1. Hemodialysis Concentrate
 - 2-1-2. Disinfectants

2-2. Dry Consumable Room


- 2-2-1. Hollow Fiber Dialyzer
- 2-2-2. Hemodialysis Blood Circuits
- 2-2-3. AVF/Dialysis Catheter Needle
- 2-2-4. Accessories Tubing for HDF
- 2-2-5. Hemodialysis Powders
- 2-2-6. Dialysis Nursing Kits
- 2-2-7. Dialysis Fluid Filter
- 2-2-8. PTA High Pressure Balloon Dilatation Catheter

3. Dialysis Furniture

- 3-1. Dialysis Bed
- 3-2. Dialysis Chair
- 3-3. Medical Trolley

4. Other Consumables

- 4-1. Dialysis Test Strip
- 4-2. Dialysis Weighting Scale
- 4-3. Other Accessories/Necessities







Technical Data Hemodialysis Machine W-T2008-B

Technical Data:

General data:

1380 x 380 x 400 mm (H x W x D); 66 KG; Dimensions

Weight approx.

Water supply:

Water inlet pressure Water inlet temperature 5°C-32°C Max. drain height

Concentrate supply

0 to 2 m suction height Supply pressure

Electrical data

Power supply AC220V±10%, 47-63HZ

Current consumption max. max. 10 A

Extracorporeal circuit

Arterial pressure monitoring

Display range Accuracy

-380 mmHg to +400 mmHg ±10 mmHg

Resolution

Venous pressure monitoring

Display range -180 mmHg to +600 mmHg ±10 mmHg

Accuracy Resolution

Transmembrane pressure monitoring -180 mmHg to +600 mmHg

Display range

±20 mmHg Resolution

Arterial blood pump

Blood flow range 30 to 600ml/min in 8 mm bloodline systems

Accuracy

Air bubble pump By ultrasound transmission, additional optical monitoring in venous clamp.

Heparin pump:

Delivery range Bolus function max. 0ml/h to 10ml/h 20 mL per bolus

Dialysis fluid circuit

Syringe size

Dialysis fluid flow range Selectable 100-800 mL/min

Dialysis fluid temperature

Accuracy

Acid concentration dialysis fluid
Default mixing ratio 1+34 (others possible)
Range 125 to 150 mmol/L

Dialysis fluid conductivity

Bicarbonate concentration dialysis fluid

Default mixing ratio

1+27.6 (others possible) -8 to +8 mmol/L bicarbonate

12 to 16 mS/cm (25°C)

±0.1 mS/cm

Ultrafiltration

0ml/h~4000ml/h UF rate Accuracy ±1% Allowed dialyzer UF factor

Parameter displayed UF goal, UF time, UF rate, UF volume

Blood leak detector

Sensitivity ≤0.35 mL blood/min (HCT=32)

KT/V Calculation

Disinfection and cleaning programs

Temperature/flow

30 to 60°C / 300 to 800ml/min

Hot Rinse (recirculation)

30 to 85°C / 300 to 800ml/min Temperature/flow

Hot disinfection with chemical

Temperature/flow 30 to 85°C / 300 to 800ml/min

Various program combination selectable

Technical Data-Options

Blood pressure Monitor (BPM)

Display range Systole: 40-280 mmHg Diastole: 40-280 mmHg Accuracy: 1 mmHg

Central Delivery system for acid conc.

1 m suction height with Central delivery system: 0-500 mbar

Endotoxin filter-- Dialysis fluid filter system ±0.1% of dialysate flow Balancing accuracy

Bicarbonate dry

Concentrate

Bi-cart SXG-F

江西三鑫医疗科技股份有限公司

JIANGXI SANXIN MEDTEC CO.,LTD.

地址: 江西南昌小蓝经济技术开发区富山大道999号

Add: No.999 Fushan Road, Xiaolan Economic Development

Zone 330200 Nanchang, Jiangxi P.R. China

电话 (Tel): 86-791-85950275 Fax: 86-791-85988819

邮箱 (E-mail): sales@sanxin-med.com 网址 (Web): www.sanxin-med.com





Hemodialysis Machine




DIASAFE Plus ensures ultrapure dialysate.

Decrease the danger of being infected by endotoxin.


Automatic test

Times of using can be recorded.

Remind of changing



Age Sex





Blood pressure monitor

Automatic NIBP monitor

Blood pressure diagram.

Measuring time setting

History record

Features:

- Intelligent Double Operation System
- LCD touch screen with button interface
- Emergency power 30Mins (Optional)
- Blood Pump&Spare Pump (for standby and also can be used for hemoperfustion).
- Heparin Pump
- Hydraulic compartment(Balance chamber + UF pump).
- Operation, Alarm information Memory function.
- A/B Ceramic proportion pump, High precision, Corrosion-proof, Accuracy.

Central Delivery System

Features:

- Nitrogen generator to inhibit bacterial growth effectively;
- Multiple safety protection: ion concentration monitoring, endotoxin filter to ensure dialysis safety;
- Filter: Filter the undissolved particles in the dialysate to achieve the hemodialysis requirement of the dialysate, effectively ensuring the quality of the concentrate.



Ro Water Purification System

Features:

- The first grade three reverse osmosis water treatment equipment in China;
- The water is more pure, the dialysis quality is higher, and the patient experience is more comfortable.



Dialyzer Reprocessing Machine



Automactic Mixing System

