

Sansin



DIALYSIS FLUID FILTER

*Pictures for reference only. In kind prevail.

- With patented technology, the specially made hollow fiber membrane has excellent biocompatibility, hydrophilicity, ultra-high flux and strong bacterial endotoxin retention capacity.
- Compatible to multi modes of dialysis machine in the market and meet diversified clinical needs.



Technical Data	
Filtering rate	Filtering rate $\geq 1000\text{mL/min}$
Particle filtration performance	Particle $\geq 10\mu\text{m/mL}$ shall less than 25pcs, Particles $\geq 25\mu\text{m/mL}$ shall less than 3pcs
Microbial filtration performance	Microbial amount $\leq 1\text{CFU/10mL}$
Endotoxin filtration performance	Endotoxin content $< 0.03\text{EU/mL}$
Shelf life	At least 150 treatments or 900 hours

Specification					
Model	Connectors		Protective Cap		Compatible Machine
	4-Con	3-Con	Press Type	Seal Ring Type	
A - I	/	✓	✓	/	Nipro
A - II	/	✓	/	✓	Nipro
A - V (01)	✓	/	✓	/	BBraun/Bellco/Toray/Wesley
A - V (02)	✓	/	✓	✓	Fresenius
A - VI	✓	/	/	✓	BBraun/Bellco/Toray/Wesley

江西三鑫医疗科技股份有限公司
JIANGXI SANXIN MEDTEC CO.,LTD.

地址: 江西南昌小蓝经济技术开发区富山大道999号
Add: No.999 Fushan Road, Xiaolan Economic Development
Zone 330200 Nanchang, Jiangxi P.R. China

电话 (Tel): 86-791-85950275 Fax: 86-791-85988819

邮箱 (E-mail): sales@sanxin-med.com

网址 (Web): www.sanxin-med.com

